

# DaVita MCoE

UI/UX Review – Falcon Physician iPad (Live Production v44.1)

January 22, 2016



## Contents

- 3 Global Findings
- 4 UI/UX Review of Falcon Physician Universal iPad (Live Production v44.1)


## Global Findings

1. **Screens are slow to load** – loading indicator appears often and slows down interaction. Could be a major usability burden when used by physicians in a fast-pasted environment **(Must Fix)**
2. **Portrait mode not supported** – users expect iPad apps to function in both portrait & landscape orientations **(Must Fix)**
3. **Inconsistent use of font styles & sizes** – should have rules and be consistent throughout. **(Highly Recommended)**
4. **Accessibility of text/button labels** – contrast ratios of text in many areas are below the minimum required contrast ratio. **(Must Fix)**



# UI/UX Review Sign In Screen


Contrast is 1.7:1  
(min is 6:1 or  
match Apple at  
3.3:1)  
**(Must Fix)**



Contrast is 3.4:1  
(min is 6:1)  
**(Must Fix)**

Logo appears to be  
low-resolution  
**(Must Fix)**

# UI/UX Review Landing Screen



Only shows recent patients by default; relies on searching to find all patients; suggest using search as a shortcut only. **(Recommended)**

Contrast is 3.3:1 (min is 6:1) **(Must Fix)**

Navigation tab is disabled; instead, should always be active and have a designed "empty state." **(Must Fix)**

Contrast is 1.3:1 (min is 6:1) **(Must Fix)**

UI/UX Review


# Patient Landing

Contrast is 3.3:1  
(min is 6:1)  
**(Must Fix)**

Extra space between  
"Blood" and  
"Pressure" **(Must Fix)**

This screen also includes Target  
Weight History, but it's not  
included in this label **(Must Fix)**

Suggest Patient Name go here  
instead of the left to make  
better use of header area  
**(Highly Recommended)**



UI/UX Review

### Patient Details: Demographics


Contrast is 3.3:1  
(min is 6:1)  
**(Must Fix)**

Contrast is 3.1:1  
(min is 6:1)  
**(Must Fix)**

Contrast is 2.5:1  
(min is 6:1 or  
match Apple at  
3.3:1)  
**(Must Fix)**

For all Patient Details screens,  
can name of data set go here  
(i.e. "Demographics") instead  
of as part of the table below?  
Perhaps re-think how both  
Patient Name & specific Detail  
view work together to make  
it as scannable, or as easily  
decipherable, as possible.  
**(Highly Recommended)**


The Patient Details views  
are paginated, but there's  
no indication. Suggest  
incorporating pagination  
element.  
**(Highly Recommended)**



UI/UX Review

Patient Details: Access

See note earlier  
in document



This interaction (refresh  
button) is a bit out-of-  
date. Suggest using  
pull-to-refresh instead  
**(Highly Recommended)**


Contrast is 2.6:1  
(min is 6:1)  
**(Must Fix)**



UI/UX Review

### Patient Details: Allergies

See note earlier in document



See note earlier in document


UI/UX Review

### Patient Details: Allergy Details

Contrast is 2.5:1  
(min is 6:1 or  
match Apple at  
3.3:1)  
**(Must Fix)**

Can this say "Allergies" to  
indicate this back button  
returns to Allergy list?  
**(Recommended)**

Can the name of the  
allergy reside here  
instead of in list?  
**(Recommended)**



UI/UX Review

Patient Details: Meds


See note earlier in document

See note earlier in document

See note earlier in document

See note earlier in document


See note earlier in document



UI/UX Review

### Patient Details: Problems

See note earlier in document



See note earlier in document


UI/UX Review

Patient Details: Problem Details

Can this say "Problems" to indicate this back button returns to Allergy list? (Recommended)


Can the name of the problem reside here instead of in list? (Recommended)

See note earlier in document



UI/UX Review


**Patient Details: Labs (refresh bug)**



Appears to be a bug here – upon first loading this screen, it appeared like this. When I tapped the Refresh button, it then seemed to load correctly. The bug did not recur the second time I navigated to this screen.  
**(Must Fix)**

UI/UX Review

Patient Details: Labs



See note earlier in document

Text size & styling different here than elsewhere in the app (Highly Recommended)

See note earlier in document

Would suggest re-thinking the design for Labs Data, i.e. shading, styling, expand/collapse interaction, size of touch targets, dynamic width of cells, etc. Also, if table slides from right to left for more data points, this could interfere with the swipe interaction to move between the Patient Details sections. (Highly Recommended)

UI/UX Review

**Patient Details:** Blood Pressure, Fluid Management, Target Weight History

See note earlier in document

See note earlier in document

See note earlier in document

Strange lines and other artifacts (Must Fix)





UI/UX Review

# Patient Details: Cultures

See note earlier  
in document



See note earlier  
in document

See note earlier  
in document

# MCoE Contacts

**Courtney Brillhart**

[courtney.brillhart@davita.com](mailto:courtney.brillhart@davita.com)

